

drop test rig plc controlled



drop test rig, plc controlled, weights from 9,1kg to 136kg, with different weight pads and seat impactor

The plc controlling allows the configuration of the test in accordance to the selected test standard. The test runs automatically with the given parameters. The selected parameters and the current condition of the test is shown on the display. The RS-485 interface can deliver data on the PC at the end of the test.

consisting of:

- plc-controlling with RS 485-interface for data transmission
- 1 test frame conceived for 136 kg
- 2 electrical stroke systems for the height adjustment of the cross traverse
- 1 drop tester stroke controlled (136 kg and 500 mm)
- 1 drop weight pad, basis weight 36 kg acc. to BIFMA
- 1 drop weight pad for arm rest-test acc. to BIFMA
- 1 additional set of weights for drop weight pads acc. to BIFMA
(bags filled with steelballs diameter 0,5 mm)
- 1 Seatimpector for drop tester with impact pad.

technical data

| | |
|--------------------------|--|
| field of application | automatic testing acc. to standards DIN 4551, ISO 7173, NEN 1812, BIFMA X5.1 and BSI/BS 5459 in development and quality control during batch production |
| test objects | chairs, revolving chairs, armchairs, stools, rocking chairs.... |
| system configuration | impact test rig for seating furniture with drop weights of 9.1 to 136 kg, seat impactor, drop heights up to 500 mm, optional: depth suspension measurement to acc. EN 581-2 |
| parameter limits | weights of 9.1 to 136 kg (depends on drop weights and additional loads) part nr. 40-910-056 drop height: up to 320 mm part nr. 40-910-057 drop height: up to 500mm number of cycles: 1 – 999 999 single cycle times: 1 – 600 sec |
| test procedure | tested pieces are subjected either to the entire test procedure prescribed by the standards or to selected individual tests. |
| supply of compressed air | air pressure: 6 bar consumption: up to 200 l/min quality of compressed air: acc. to ISO/DIS 8573-1 - residual oil content better class 1 - residual dust content better class 1 - residual water content better class 4 |
| input control unit | 230 V/50 Hz |
| dimensions of test rig | 1500 x 1700 x 2850 (length x width x height) |
| weight | 450 kg |
| optional accessories | measurement device for the finding of the acceleration of the drop weight acc. to EN 581-2 |

test equipment
furniture test rigs
type 40-910-056/057



Hegewald & Peschke
Meß- und Prüftechnik GmbH

Hegewald & Peschke
Meß- und Prüftechnik GmbH
Am Gründchen 1 • 01683 Nossen
Telefon: +49 (0) 3 52 42 - 4 45 10
Telefax: +49 (0) 3 52 42 - 4 45 11
mail: info@hegewald-peschke.com
<http://www.hegewald-peschke.com>



Your partner:
